



Vivaspin 6 and 20 ml

Technical data and operating instructions.
For in vitro use only.

New

3 kDa MWCO
in Vivaspin 6

Vivascience Vivaspin 6 and 20 ml

Storage conditions / shelf life

Vivaspin ultrafiltration spin columns should be stored at room temperature.

The devices should be used before the expiry date printed on the box.

Introduction

Vivaspin concentrators are disposable ultrafiltration devices for the concentration and/or purification of biological samples. Vivaspin 6 is suitable for sample volumes of 2-6 ml and the Vivaspin 20 can handle samples up to 20 ml. Both products feature twin vertical membranes for unparalleled speed.

Vivaspin 20 purification alternatives include a diafiltration cup that allows one step removal of salts and other contaminating micromolecules, and a gas pressure mode for increased flexibility and even faster processing.

The innovative design (US Patent No. 5,647,990, second patent pending), ease of use, speed and exceptional concentrate recoveries are the main features of the concentrators.

Centrifugal Operation

Vivaspin concentrators can be used in swing bucket or fixed angle rotors accepting standard conical bottom tubes. In a single spin, solutions can be concentrated in excess of 100 X. Samples are typically concentrated in 10 to 30 minutes with macromolecular recoveries in excess of 95 %.



The longitudinal membrane orientation and thin channel concentration chamber, provide optimum cross flow conditions even for particle laden solutions; the centrifugal force pulling particles and solids away from the membrane to the bottom of the device. Macromolecules collect in an impermeable concentrate pocket integrally moulded below the membrane surface, thereby eliminating the risk of filtration to dryness.

Pressurised Operation

When an appropriate centrifuge is unavailable, or for single sample processing, Vivaspin 20 can be filled with up to 15 ml and pressurised for bench top concentration. For even faster processing, pressure can be combined with centrifugal force. "Pressure-Fugation" is particularly suitable for viscous samples such as serum, or when processing at low temperatures, and generally when minimum process time is essential.

Equipment Required

A. For use with centrifuge

1. Centrifuge with swing bucket or fixed angle rotor (minimum 25°).
2. Pasteur or fixed volume pipettes for sample delivery and removal.

| Device | Carrier Required |
|-------------|------------------|
| Vivaspin 6 | 15 ml / 17 mm Ø |
| Vivaspin 20 | 50 ml / 30 mm Ø |

B. For use with Pressure (Vivaspin 20 only)

1. Vivaspin 20 Pressure Head (Product No. VCA200).
2. Charge Valve for Pressure Head (Product No. VCA005).
3. Vivascience Air Pressure Controller (Product No. VCA002). or equivalent pressure regulator

For use with Pressure and Centrifuge

1. All of the equipment shown in A. and B. above.

| Equipment required | Vivaspin 6 | | Vivaspin 20 | |
|---|---|-------------|---|-------------|
| Centrifuge | | | | |
| Rotor type | Swing bucket | Fixed angle | Swing bucket | Fixed angle |
| Minimum rotor angle | - | 25° | - | 25° |
| Rotor cavity | To fit 15 ml (17 mm) conical bottom tubes | | To fit 50 ml (30 mm) conical bottom tubes | |
| Optional pressure accessories for Vivaspin 20 | | | | |
| Air pressure controller (APC) complete with pressure gauge, regulator, over-pressure safety valve, female connector to Vivascience pressure products and 1m extension line (4 mm pneumatic tubing) with male and female connectors and 1 m of 6 mm inlet tubing | | | Prod no. VCA002 | |
| Charge valve | | | Prod no. VCA005 | |
| VS20 pressure head | | | Prod. no.VCA200 | |
| Concentrate recovery | | | | |
| Pipette type | Fixed or variable volume | | Fixed or variable volume | |
| Recommended tip | Thin gel loader type | | Thin gel loader type | |

Rotor compatibility

Please note: Vivaspin 20 (30 mm x 116 mm) is designed to fit into rotors that can accommodate Falcon 50 ml conical bottom tubes, e.g. Beckman Allegra 25R with TS-5.1-500 swing-out rotor with BUC 5 buckets and 368327 adaptors; Beckman TA-10.250 25° fixed angle rotor with 356966 adaptors; Heraeus Multifuge 3 S-R with (Heraeus / Sorvall) 75006445 swing out rotor with 75006441 buckets and adaptors for Falcon 50 ml conical bottom tubes.

These devices are not designed to fit into rotors that only accept round bottom 29 mm x 105 mm tubes, e.g. Sorvall SS34 or Beckmann JA 20. If your rotor accepts only 29 mm x 105 mm round bottom tubes, please use the Vivaspin 15, which can be used in either round bottom or conical centrifuge tubes.



Operation

In Centrifuge

1. Select the most appropriate membrane cut-off for your sample.
For maximum recovery select a MWCO at least 50 % smaller than the molecular size of the species of interest.
2. Fill concentrator with up to maximum volumes shown in table 1.
(Ensure screw closure is fully seated)
3. Insert assembled concentrator into centrifuge (when fixed angle rotors are used, angle concentrator so that the printed window faces upwards/outwards).
4. Centrifuge at speeds recommended in table 2, taking care not to exceed the maximum g force indicated by membrane type and MWCO
5. Once the desired concentration is achieved, (see table 3a and 3b for guide to concentration times), remove assembly and recover sample from the bottom of the concentrate pocket with a pipette.

Using Gas Pressure (Vivaspin 20 only)

1. Select appropriate membrane as above.
2. Fill concentrator (maximum 15 ml).
3. Fit Pressure Head (Prod. No. VCA200), and hand tighten to ensure an air tight seal.
4. Using Air Pressure Controller (APC), exchange the Female Coupling for the Charge Valve (Prod. No.VCA005), on the APC Extension line. Pressurise by pressing the Charge Valve into the inlet valve of the pressure head to achieve an air tight seal.
5. EITHER - stand assembly in a rack and allow to concentrate,
OR - for faster processing, insert assembled, pressurised concentrator into centrifuge and spin (see table 1 for spin speeds).
6. Once the desired concentration has been reached, (see table 3b for guide to concentration times), remove assembly and de-pressurise by unscrewing cap.
7. Extract concentrate using a pasteur or fixed volume pipette.



Desalting/Buffer Exchange

1. Concentrate sample to desired level.
2. Empty filtrate container.
3. Refill concentrator with an appropriate solvent.
4. Concentrate the sample again and repeat the process until the concentration of contaminating microsolite is sufficiently reduced. Typically, 3 wash cycles will remove 99 % of initial salt content.

One Step Desalting with Vivaspin 20

Salts and contaminants can be removed in a single step when using the special diafiltration cup available with the Vivaspin 20. This is due to the constant washing action (constant volume dialfiltration), of the buffer solution in the cup as it replaces solvent and salts passing through the ultrafiltration membrane.

1. Place 2nd sample solution in the concentrator. (Larger volumes can be desalted by first concentrating down to 2 ml and decanting filtrate).
2. Empty filtrate container.
3. Insert diafiltration cup into concentrator and fill with 10ml deionised water or buffer solution. Re-fit blue lid over the diafiltration cup.
4. Repeat concentration process; over 98 % of salts will be removed in this step.
5. Remove diafiltration cup and recover concentrated and purified sample.

Usage Tips

1. Flow Rate

Filtration rate is affected by several parameters, including MWCO, porosity, sample concentration, viscosity, centrifugal force and temperature. Expect significantly longer spin times for starting solutions with over 5 % solids. When operating at 4° C. flow rates are approximately 1.5 times slower than at 25° C. Viscous solutions such as 50 % glycerine will take up to 5 times longer to concentrate than samples in a predominantly buffer solution.

2. Prerinsing

Membranes fitted to Vivaspin concentrators contain trace amounts of Glycerine and Sodium azide. Should these interfere with analysis they can be removed by rinsing fill volume of buffer solution or deionised water through the concentrator. Decant filtrate and concentrate before processing sample solution.

3. Sterilisation of

Polyethersulfone Membranes

Polyethersulfone membranes should not be autoclaved as high temperatures will substantially increase membrane MWCO. To sterilise, use a 70 % ethanol solution or sterilising gas mixture.

4. Chemical Compatibility

Vivaspin concentrators are designed for use with biological fluids and aqueous solutions. For chemical compatibility details, refer to table 4.

| Table 1: Technical specifications | Vivaspin 6 | Vivaspin 20 |
|-----------------------------------|---------------------|---------------------------|
| Concentrator capacity | | |
| Swing bucket rotor | 6 ml | 20 ml |
| Fixed angle rotor | 6 ml | 14 ml |
| With pressure head | - | 15 ml |
| Dimensions | | |
| Total length | 122 mm | 116 mm |
| | - | 125 mm with pressure head |
| Width | 17 mm | 30 mm |
| Active membrane area | 2.5 cm ² | 6.0 cm ² |
| Hold up volume of membrane | <10 µl | < 20 µl |
| Dead stop volume | 30 µl | 50 µl |
| Materials of construction | | |
| Body | Polycarbonate | Polycarbonate |
| Filtrate vessel | Polycarbonate | Polycarbonate |
| Concentrator cap | Polypropylene | Polypropylene |
| Pressure head | - | Acetal/aluminium |
| Membrane | Polyethersulfone | Polyethersulfone |

| Table 2: Recommended Spin Speed (xg) | | | |
|--------------------------------------|--------------|---------------|--------------------------|
| Vivaspin 6 | Swing Bucket | Fixed Angle | |
| Membrane | max | max | |
| 3 - 50,000 MWCO PES | 4,000 | 10,000 | |
| >100,000 MWCO PES | 4,000 | 6,000 | |
| Vivaspin 20 | Centrifuge | Pressure-Fuge | |
| Rotor | Swing Bucket | Fixed Angle | Swing Bucket (5 bar max) |
| Membrane | max | max | max |
| 3 - 50,000 MWCO PES | 5,000 | 8,000 | 3,000 |
| >100 - 300,000 MWCO PES | 3,000 | 6,000 | 2,000 |

| Table 3a: Typical performance Vivaspin 6 | | Time to concentrate 30x min. at 20°C and solute recovery % | | | |
|--|----------|--|------|-----------------|------|
| Rotor | | Swing bucket | | 25° Fixed angle | |
| Start volume | | 6 ml | | 6 ml | |
| | | Min. | Rec. | Min. | Rec. |
| Cytochrome c 0.25 mg/ml (12,400 MW) | | | | | |
| 3,000 | MWCO PES | - | - | 90 | 97% |
| BSA 1.0 mg/ml (66,000 MW) | | | | | |
| 5,000 | MWCO PES | 20 | 98% | 12 | 98% |
| 10,000 | MWCO PES | 13 | 98% | 10 | 98% |
| 30,000 | MWCO PES | 12 | 98% | 9 | 97% |
| IgG 0.25 mg/ml (160,000 MW) | | | | | |
| 30,000 | MWCO PES | 18 | 96% | 15 | 95% |
| 50,000 | MWCO PES | 17 | 96% | 14 | 95% |
| 100,000 | MWCO PES | 15 | 91% | 12 | 91% |
| Latex beads 0.004 % in DMEM +10 % FCS (0.055 µm) | | | | | |
| 300,000 | MWCO PES | - | - | 25 | 99% |
| Latex beads 0.004 % in DMEM +10 % FCS (0.24 µm) | | | | | |
| 1,000,000 | MWCO PES | - | - | 4 | 99% |
| Yeast 1.0 mg/ml (<i>S. Cerevisiae</i>) | | | | | |
| 0.2 µm | PES | 4 | 97 % | 3 | 97 % |

| Table 3b: Typical performance Vivaspin 20 | | Time to concentrate 30x min. at 20°C and solute recovery % | | | | | | | |
|--|----------|--|------|-----------------|------|-----------|------|--------------|------|
| Mode | | Centrifuge | | Centrifuge | | Bench top | | Press-fuge | |
| Rotor | | Swing bucket | | 25° Fixed angle | | Pressure | | Swing bucket | |
| Start volume | | 20 ml | | 14 ml | | 10 ml | | 10 ml | |
| | | Min. | Rec. | Min. | Rec. | Min. | Rec. | Min. | Rec. |
| Cytochrome c 0.25 mg/ml (12,400 MW) | | | | | | | | | |
| 3,000 | MWCO PES | 110 | 97 % | 180 | 96 % | 60 | 96 % | - | - |
| BSA 1.0 mg/ml (66,000 MW) | | | | | | | | | |
| 5,000 | MWCO PES | 23 | 99 % | 29 | 99 % | 50 | 98 % | 14 | 98 % |
| 10,000 | MWCO PES | 16 | 98 % | 17 | 98 % | 32 | 97 % | 8 | 97 % |
| 30,000 | MWCO PES | 13 | 98 % | 15 | 98 % | 32 | 97 % | 8 | 97 % |
| IgG 0.25 mg/ml (160,000 MW) | | | | | | | | | |
| 30,000 | MWCO PES | 27 | 97 % | 20 | 95 % | 46 | 94 % | 13 | 97 % |
| 50,000 | MWCO PES | 27 | 96 % | 22 | 95 % | 46 | 93 % | 13 | 96 % |
| 100,000 | MWCO PES | 25 | 91 % | 20 | 90 % | 42 | 88 % | 12 | 94 % |
| Latex beads 0.004 % in DMEM +10 % FCS (0.055 µm) | | | | | | | | | |
| 300,000 | MWCO PES | 20 | 99 % | 35 | 99 % | 10 | 99 % | - | - |
| Latex beads 0.004 % in DMEM +10 % FCS (0.24 µm) | | | | | | | | | |
| 1,000,000 | MWCO PES | 4 | 99 % | 12 | 99 % | 4 | 99 % | - | - |
| Yeast 1.0 mg/ml (<i>S. Cerevisiae</i>) | | | | | | | | | |
| 0.2 µm | PES | 15 | 95 % | 5 | 95 % | 20 | 95 % | 2 | 95 % |

Table 4: Chemical Compatibility (2hr contact time)

| Solutions - Compatible pH range pH 1-9 | | Solutions - Compatible pH range pH 1-9 | |
|--|----|--|----|
| Acetic Acid (25.0 %) | OK | Lactic Acid (5.0 %) | OK |
| Acetone (10.0 %) | NO | Mercaptoethanol (10 mM) | OK |
| Acetonitrile (10.0 %) | NO | Methanol (60 %) | ? |
| Ammonium Hydroxide (5.0 %) | ? | Nitric Acid (10.0 %) | OK |
| Ammonium Sulphate (saturated) | OK | Phenol (1.0 %) | ? |
| Benzene (100 %) | NO | Phosphate Buffer (1.0 M) | OK |
| n - Butanol (70 %) | OK | Polyethylene Glycol (10 %) | OK |
| Chloroform (1.0 %) | NO | Pyridine (100 %) | ? |
| Dimethyl Formamide (10.0 %) | ? | Sodium Carbonate (20 %) | ? |
| Dimethyl Sulfoxide (5.0 %) | OK | Sodium Deoxycholate (5.0 %) | OK |
| Ethanol (70.0 %) | OK | Sodium Dodecylsulfate (0.1 M) | OK |
| Ethyl Acetate (100 %) | NO | Sodium Hydroxide | NO |
| Formaldehyde (30 %) | OK | Sodium Hypochlorite (200 ppm) | ? |
| Formic Acid (5.0 %) | OK | Sodium Nitrate (1.0 %) | OK |
| Glycerine (70 %) | OK | Sulfamic Acid (5.0 %) | OK |
| Guanidine HCl (6 M) | OK | Tetrahydrofuran (5.0 %) | NO |
| Hydrocarbons, aromatic | NO | Toluene (1.0 %) | NO |
| Hydrocarbons, chlorinated | NO | Trifluoroacetic Acid (10 %) | OK |
| Hydrochloric Acid (1 M) | OK | Tween 20 (0.1 %) | OK |
| Imidazole (500 mM) | OK | Triton X-100 (0.1 %) | OK |
| Isopropanol (70 %) | OK | Urea (8 M) | OK |

OK = Acceptable ? = Questionable NO = Not recommended

Ordering information

| Vivaspin 6 Polyethersulfone | Pack size | Prod. no. | Vivaspin 20 Polyethersulfone | Pack size | Prod. no. |
|--|-----------|-----------|--|-----------|-----------|
| 3,000 MWCO | 25 | VS0691 | 3,000 MWCO | 12 | VS2091 |
| 3,000 MWCO | 100 | VS0692 | 3,000 MWCO | 48 | VS2092 |
| 5,000 MWCO | 25 | VS0611 | 5,000 MWCO | 12 | VS2011 |
| 5,000 MWCO | 100 | VS0612 | 5,000 MWCO | 48 | VS2012 |
| 10,000 MWCO | 25 | VS0601 | 10,000 MWCO | 12 | VS2001 |
| 10,000 MWCO | 100 | VS0602 | 10,000 MWCO | 48 | VS2002 |
| 30,000 MWCO | 25 | VS0621 | 30,000 MWCO | 12 | VS2021 |
| 30,000 MWCO | 100 | VS0622 | 30,000 MWCO | 48 | VS2022 |
| 50,000 MWCO | 25 | VS0631 | 50,000 MWCO | 12 | VS2031 |
| 50,000 MWCO | 100 | VS0632 | 50,000 MWCO | 48 | VS2032 |
| 100,000 MWCO | 25 | VS0641 | 100,000 MWCO | 12 | VS2041 |
| 100,000 MWCO | 100 | VS0642 | 100,000 MWCO | 48 | VS2042 |
| 300,000 MWCO | 25 | VS0651 | 300,000 MWCO | 12 | VS2051 |
| 300,000 MWCO | 100 | VS0652 | 300,000 MWCO | 48 | VS2052 |
| 1,000,000 MWCO | 25 | VS0661 | 1,000,000 MWCO | 12 | VS2061 |
| 1,000,000 MWCO | 100 | VS0662 | 1,000,000 MWCO | 48 | VS2062 |
| 0.2 µm | 25 | VS0671 | 0.2 µm | 12 | VS2071 |
| 0.2 µm | 100 | VS0672 | 0.2 µm | 48 | VS2072 |
| Starter pack (5 of each 5 k, 10 k, 30 k, 50 k, 100 k) | 25 | VS06S1 | Starter pack (2 of each 5 k, 10 k, 30 k, 50 k, 100 k, 0.2 µm) | 12 | VS20S1 |
| | | | Vivaspin 20 accessories | | |
| | | | Air pressure controller (APC) | 1 | VCA002 |
| | | | Charge valve for pressure head | 1 | VCA005 |
| | | | Diafiltration cups | 12 | VSA005 |
| | | | Female connector | 1 | VCA010 |
| | | | Male connector | 1 | VCA011 |
| | | | 4 mm OD pneumatic tube (3 m) | 1 | VCA012 |
| | | | Vivaspin 20 pressure head | 1 | VCA200 |

Vivascience ultrafiltration product range at a glance

| Product | Sample volume | Mode | Membranes available |
|----------------|-----------------------|-----------------------------|---|
| Vivaspin 500 | 100 µl - 600 µl | Centrifugal | Polyethersulfone |
| Vivaspin 2 | 0.4 ml - 2 ml | Centrifugal | Polyethersulfone, Cellulose Triacetate Regenerated Cellulose, Hydrosart® |
| Centrisart | 0.5 ml - 2.5 ml | Centrifugal | Polyethersulfone, Cellulose Triacetate |
| Vivaspin 4 | 1 ml - 4 ml | Centrifugal | Polyethersulfone |
| Vivaspin 6 | 2 ml - 6 ml | Centrifugal | Polyethersulfone |
| Vivaspin 15 | 2 ml - 15 ml | Centrifugal | Polyethersulfone |
| Vivaspin 15R | 2 ml - 15 ml | Centrifugal | Hydrosart® |
| Vivaspin 20 | 5 ml - 20 ml | Centrifugal Gas pressure | Polyethersulfone |
| Vivacell 70 | 10 ml - 70 ml | Centrifugal Gas pressure | Polyethersulfone |
| Vivacell 100 | 20 ml - 100 ml | Centrifugal Gas pressure | Polyethersulfone |
| Vivacell 250 | 50 ml - 250 ml | Gas pressure | Polyethersulfone |
| Vivaflow 50 | 100 ml - > 5 l | Tangential flow | Polyethersulfone, Regenerated Cellulose |
| Vivaflow 200 | 500 ml - > 5 l | Tangential flow | Polyethersulfone, Regenerated Cellulose, Hydrosart® |
| Vivapore 2 | 0.5 ml - 2.5 ml/15 ml | Solvent absorption | Polyethersulfone, Regenerated Cellulose |
| Vivapore 5 | 1 ml - 5 ml | Solvent absorption | Polyethersulfone, Regenerated Cellulose |
| Vivapore 10/20 | 2 ml - 10 ml/20 ml | Solvent absorption | Polyethersulfone, Regenerated Cellulose |
| Vivapore Q5 | 0.5 ml - 5 ml | Solvent absorption | Polyethersulfone |
| Vivapore Q10 | 1 ml - 10 ml | Solvent absorption | Polyethersulfone |



New Vivapure kits and devices now available

The development of the innovative Vivapure range continues with the launch of new protein purification kits and devices. Vivapure spin columns are currently available with metal chelate protein A, epoxy and ion exchange membrane chemistries.

Contact us for more details or visit www.vivascience.com

| Vivascience Technical Support | | Phone | Fax | E-mail |
|-------------------------------|---|-------------------------------|--------------------|------------------------------|
| USA | Vivascience Service & Technical Support | +1 877 452 2345 (toll free) | +1 631 253 5460 | info.usa@vivascience.com |
| Europe | Vivascience Support Center | +49 1 802 848 201 (toll free) | +49 1 802 848 202 | info@vivascience.com |
| International | Vivascience Support Center | +49 511 524 875 60 | +49 511 524 875 69 | info@vivascience.com |
| Vivascience Customer Sales | | | | |
| France | Vivascience S.A.R.L. | +33 169 19 93 23 | +33 160 13 95 05 | info.france@vivascience.com |
| Germany | Vivascience AG | +49 551 308 4023 | +49 551 308 3289 | info.germany@vivascience.com |
| UK | Vivascience Ltd. | +44 1372 737 159 | +44 1372 726 171 | info.uk@vivascience.com |
| USA | Vivascience Inc. | +1 877 452 2345 | +1 631 253 5445 | info.usa@vivascience.com |

Vivascience AG

Feodor-Lynen-Strasse 21 Phone: +49 511 524875-0 E-mail: info@vivascience.com
 30625 Hannover, Germany Fax: +49 511 524875-19 Web: www.vivascience.com

The information in this publication is believed to be accurate, and is offered in good faith for the benefit of the user. As we cannot anticipate or control conditions of product application, we do not warrant suitability or favourable results, nor that use will not infringe on any patents. Products are furnished upon the condition that the user assumes all risks and liabilities and that neither seller or manufacturer shall be liable for any loss or damage, direct or consequential arising from the use or inability to use the product. Specifications subject to change without notice. © 2002 Vivascience Ltd. Publication No. SLU6092-e05019. Order No. 85030-511-57