

RNA extraction protocols

The RNA extraction of cells cultured in BIOMIMESYS® is performed according to the manufacturer's instructions of commercial mRNA extraction kits.

BIOMIMESYS® is compatible with manual (Trizol) and commercially available mRNA extraction kits such as QIAGEN RNeasy Micro or Macherey-Nagel NucleoSpin RNA XS.

Below are recommendations to optimise RNA extraction:



- NucleoSpin RNA XS kit from Macherey-Nagel:
- Lyse cells by vortexing hydrogels thoroughly 2 x 5 seconds in 100μL buffer RA1 + 2μL TCEP (for 2 hydrogels).
- Directly place the 2 hydrogels, side by side, in the purple clarification column, then place iton a collecting tube.
- Centrifuge 5 minutes at 10 000g.
- Measure the recovered filtrate volume with a pipette and add the same volume of 70% ethanol. Mix by pipetting (at least 5 times) before dispensing the mixture onto the blue column and continue RNA extraction according to the manufacturer instructions (step 6).



- RNeasy Micro Kit from Qiagen :
- Lyse cells by vortexing hydrogels thoroughly 2 x 5 seconds in 350μL buffer RLT (for 2 hydrogels).
- Break down hydrogels in lysis buffer and by aspirating total volume in 1mL syringe with a 1.9mm needle.
- Centrifuge 5 minutes at 10 000g, collect supernatant.
- Measure the collected supernatant volume, add same volume of 70% ethanol and mix well by pipetting several times.
- Transfer the mixture onto the RNeasy Min Elute spin (pink column) and continue RNA extraction according to the manufacturer's instructions (step 3 from Quick Start Protocol).



V2-12 201^e

BIOMIMESYS*

RNA extraction protocols

analytikjena • InnuPrep RNA Mini from AnalytikJena :

- Lyse cells by vortexing hydrogels thoroughly 2 x 5 seconds in 400μL buffer RL (for 2 hydrogels).
- Directly place the 2 hydrogels, side by side, in the blue Spin Filter D column, and place it on a collecting tube.
- Centrifuge 5 minutes at 10 000g.
- Measure the filtrate volume, add same volume of 70% ethanol and mix well by pipetting several times.
- > Transfer the mix onto the Spin Filter R column (purple column) and continue RNA extraction according to the manufacturer's instructions (step 4 RNA extraction from eucaryotic cells protocol).
- Using Trizol/Chloroform/Isopropanol extraction :
 - Lyse cells by vortexing hydrogels thoroughly 2 x 5 seconds in 1ml Trizol (for 1 to 5 hydrogels).
 - Break down/destroy hydrogels in Trizol by aspirating total volume in 1mL syringe with 1.9mm needle.
 - Optional step: successive freeze-drying of lysates containing hydrogels.
 - Add 200μL of Chloroform, vortex 2 x 5 seconds and incubate 3 minutes at room temperature.
 - Centrifuge 20 minutes at 10 000g at 4°C.
 - Collect aqueous phase in a new tube.
 - > Add 500μL of Isopropanol, vortex 2 x 5 seconds and incubate 20 minutes at room temperature or overnight at -20°C.
 - Centrifuge 20 minutes at 10 000g at 4°C.
 - Discard supernatant without aspirating the pellet, add 1 mL of 75% ethanol and vortex
 2 x 5 seconds.
 - Centrifuge 15 minutes at 10 000g at 4°C.
 - Remove supernatant without aspirating the pellet, add 1 mL of 75% ethanol, vortex 2 x 5 seconds, and centrifuge 15 minutes at 10 000g at 4°C.



V2-12 201^e



RNA extraction protocols

- Vacuum or air dry the RNA pellet.
- Resuspend the RNA pellet in 10-20μL RNase-free water by passing the solution up and down several times through a pipette tip.
- ➤ Incubate in a water bath or heat block set at 55°C for 5 minutes.
- > Purify obtained RNA or proceed to downstream applications.

Yield per hydrogel of a 96 well plate, according to kits used:

Kits	Cells	Culture time	Yield/hydrogel	Ratio (OD 260/280nm)
NucleoSpin RNA XS Macherey- Nagel	HT-29	7 days	1-2 µg	2.09-2.12
		14 days	1.5-3 μg	
		30 days	8 µg	
	HepG2	30 days	600 ng	2.09-2.08
	3T3-L1	8 days	840 ng	2.14-2.18
innuPrep RNA Mini Eurobio	HT-29	30 days	8 µg	2.09-2.13
	HepG2	30 days	400 ng	
RNeasy Micro Kit Qiagen	3T3-L1	8 days	940 ng	2.08-2.09
Trizol/ Chloroform/ Isopropanol	HT-29	30 days	14-23 µg	1.83-1.94
	HWP	6 days	681 ng	1,83-2,13
		12 days	411 ng	
		19 days	417 ng	
		26 days	593 ng	
		31 days	491 ng	



V2-12.201