



Protein Quantification: Bradford Protein Reagent from Fermentas

Ready-to-use formulation: no dilution is required

Rapid colour development

Broad working concentration range: 1-2,000 µg/ml

Compatible with a variety of substances present in samples, including solvents, salts, detergents, thiols, reducing and chelating agents (full chemical compatibility can be found on www.generon.co.uk)



Bradford 1 x Reagent-Ready to use

Cat no: **GEN-R1271**

Rapid and accurate estimation of protein concentration in solution

Assay can be performed in 1.5/2.2 ml microfuge tubes or multiwell plates

Standard curve obtained using Bradford reagent, ready-to-use, in a standard assay (Fig. 1) and in a microassay (Fig. 2)

Enough reagent for 8,000 microassays in a multiwell plate

BSA - bovine serum albumin
BGG - bovine gamma globulin

Ordering Information:

Product	Amount	Code	Price
Bradford Reagent	1 L	GEN-R1271	£57.00



generon

Tel: (01753) 866-511

Fax: (01753) 864-942

E-mail: info@generon.co.uk

Generon Ltd,
12 Rawcliffe House, Howarth Road,
Maidenhead,
Berkshire SL6 1AP, UK.
Web: www.generon.co.uk