

## Electrophoresis and Western Blots: Protein Ladders from Fermentas


- Sharp bands
- Stable colours
- Broad range
- Multicoloured bands for easy orientation
- Easy-to-remember pattern
- Ready-to-use and supplied in a loading buffer for direct loading
- Stable at +4°C for 3 months;  
Long term storage at -20°C

PageRuler™ Prestained Protein Ladder

PageRuler™ Plus Prestained Protein Ladder

Spectra™ Multicolour Broad Range Protein Ladder

PageRuler™ Unstained Protein Ladder



### PageRuler™ Prestained Protein Ladder

Cat no: **GEN-SM0671**


10 proteins in total

Green 10 kDa protein for easy orientation

10-170 kDa protein ladder

2 x 250 µl volume

50-100 applications (5 µl/well)



### PageRuler™ Plus Prestained Protein Ladder

Cat no: **GEN-SM1811**


9 proteins in total

Two orange proteins for easy orientation

10-250 kDa protein ladder

2 x 250 µl volume

50-100 applications (5 µl/well)



### Spectra™ Prestained Protein Ladder

Cat no: **GEN-SM1841**


10 proteins in total

Multicoloured proteins for easy orientation

10-260 kDa protein ladder

2 x 250 µl volume

25-50 applications (10 µl/well)



### PageRuler™ Unstained Protein Ladder

Cat no: **GEN-SM0661**

14 proteins in total

10-200 kDa protein ladder

2 x 250 µl volume

50-100 applications (5 µl/well)

### Ordering Information:

Product	Range	Amount	Code	Price
PageRuler™ Prestained Ladder	10-170 kDa	2 x 250 µl	GEN-SM0671	£80.50
PageRuler™ Plus Prestained Ladder	10-250 kDa	2 x 250 µl	GEN-SM1811	£87.50
Spectra™ Prestained Ladder	10-260 kDa	2 x 250 µl	GEN-SM1841	£98.50
PageRuler™ Unstained Ladder	10-200 kDa	2 x 250 µl	GEN-SM0661	£64.50

# generon

Tel: (01753) 866-511

Fax: (01753) 864-942

E-mail: [info@generon.co.uk](mailto:info@generon.co.uk)

**Generon Ltd,**  
**12 Rawcliffe House, Howarth Road,**  
**Maidenhead,**  
**Berkshire SL6 1AP, UK.**  
**Web: [www.generon.co.uk](http://www.generon.co.uk)**