



# generon

## Thompson Single StEP chromatography columns



Thompson's Single StEP Empty Columns allow you the opportunity to fill different sized columns with resins for protein purification applications. The sample is mixed with the resin and batch incubated. The columns can then be attached to an FPLC or peristaltic pump and the purified protein can then eluted with a gradient.

### Specifications:

Sinter type:	Ultra High Density Polyethylene
Plastic construction:	Polypropylene
Pore size:	25 $\mu$ m
Mode of action:	Gravity or Batch Incubation/Pump Elution
Reservoir volume:	10 ml, 25 ml, 50 ml, 100 ml, 200 ml, 320 ml, 650 ml
Storage:	Store at RT (non-sterile)

Order Code	Description	Pack size	Price
9452086-10	10 ml Single StEP column with bottom frit	10	£88.00
9452088-10	25 ml Single StEP column with bottom frit	10	£88.00
9452090-10	50 ml Single StEP column with bottom frit	10	£98.00
9452092-10	100ml Single StEP column with bottom frit	10	£130.00
9452094-10	200 ml Single StEP column with bottom frit	10	£130.00
9452099-10	320 ml Single StEP column with bottom frit	10	£300.00
9452097-10	650 ml Single StEP column with bottom frit	10	£345.00
235008-20	Dual ended Caps (Blue)	20	£20.00

# generon

Tel: (01753) 866-511  
Fax: (01753) 20-88-99  
E-mail: [info@generon.co.uk](mailto:info@generon.co.uk)

**Generon Ltd,**  
**12 Rawcliffe House, Howarth Road,**  
**Maidenhead,**  
**Berkshire SL6 1AP, UK.**  
**Web: [www.generon.co.uk](http://www.generon.co.uk)**